

Table 4. MS parameters in PRM mode of the five target analytes and their detectability in urine and plasma samples

Analyte	Precursor Ion, <i>m/z</i>	Product Ion, <i>m/z</i>	NCE, eV	RT (min)	Detection urine	Detection plasma
Unknown	220.1543	202.1436	10	2.25	Yes	No
Unknown	416.3159	98.0967	65	3.90	Yes	Yes*
Unknown	432.3108	98.0967	65	3.84	Yes*	No
Unknown	444.3108	98.0967	65	4.09	Yes*	Yes*
Unknown	597.3382	84.0812	50	4.96	Yes	Yes

*Except in PM samples